博昊科技有限公司
BOUGH TECHNOLOGY LIMITED

Unit513,Lakeside 1,Bldg.15,No.8 Science Park West Avenue,
HongKong Science Park,New Territories,HongKong.
Tel:+852 3595 1823. Fax:+852 2661 9845. www.bough.com.cn

BG544

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Ver:1.0

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Summary

This product adopts EMS technology (the working voltage can reach 30V - 60V), and the product is not used by patients with heart disease, hypertension and other diseases. Please keep this product away from the human brain. Please do not provide this product to children under 18 years of age and children. This product mainly applies to the family.

Method of application

Download APP: products provide ways to download the corresponding APP automatically connected bluetooth: open the bluetooth and open the APP, open the APP main interface, fitness apparatus after startup will itself with mobile phone bluetooth connection (need not enter the bluetooth connection password) boot: long press + boot (about 3 s) shutdown: long press - button to turn it off (about 3 s) open the APP, open the bluetooth connection abdominal massage device, and set the gear strength and work mode gear strength to adjust: fitness apparatus after startup, can press through short on a host of + and - keys to adjust the gear strength, strength of also can adjust the gear through the APP working mode switching: through long press on the host + 3 s to switch working pattern, also can switch work mode through the APP. Charge: use USB charging line to abdominal massage instrument replacement gel stick: gel belongs to consumable, use for a long time can lead to gel aging damage, influence experience effect, for the sake of safety and experience effect, please use the company specializing in the gel stick for replacement using a reset button: when fitness instrument crashes or unable to boot, the user can be reset by reset button will abdominal health meter to check the record

View log

You can view the use record in the fitness meter APP.

Solutions to anomalies

Q:Unable to connect to the APP

A: fitness reasons:

1. disconnect bluetooth, open again, make sure is connected
2. If it is Android mobile phone, intelligent system version is not lower than 4.3
3. fitness apparatus and mobile distance too far
4. Obstacles such as A wall between equipment can't normal control

A:solution:

1. open the bluetooth
2. version of the mobile phone use over A Android4.3 (Android version 6.0 and above phone need to open the GPS)
3. fitness apparatus near the phone
4. Avoid using battery in the case of there are obstacles, charging starts again after

Q: unable to start up

A: reasons:

1. lack of battery

-
2. Crash:
A: solution
1. Will be charge by fitness apparatus
 2. Will be reset by fitness apparatus

Cleaning and maintenance and storage

This machine if there is dirt, please use soft dry cloth to wipe the pollution is very serious, please use cloth dipped in water or neutral detergent and twist dry wipe, then use dry cloth to wipe clean, prevent water from entering the product, so as not to damage the internal electronic components

- After using gel belongs to consumable, please ensure cleanliness and protection
- Do not use this product exposure to strong sunlight, dust and pollution of the environment
- Please do not use gasoline, naphtha, thinner, alcohol and other chemical solvents to clean this product

Matters need attention

1. This product is a precision instrument, please handle with care, avoid collisions, fall and cause product failure
2. Do not use this product into the water
3. Such as product can't normal use, please stop using and consult the dealer
4. Please do not remove this product without authorization
5. Provided all the pictures in this manual are for reference only, the kind prevail